



14-3-3 epsilon Rabbit pAb
Catalog NO.: EA154
For research use only.

Overview

Product name	14-3-3 epsilon Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species reactivity	Rat, Mouse
Recommended dilutions	WesternBlot:1/1000 Immunohistochemistry:1/200-500 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Rabbit anti-human 14-3-3 epsilon pAb
Species	Human
Storage	PBS, 0.02% sodium azide, 50% glycerol, pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	28kDa
GeneID Human	7531
Human Swiss-Prot No.	P62258
Cellular localization	Cytoplasm, Nucleus
Alternatin	



W b a a 1) M B a , 2) Ra B a , 14-3-3
Rabb Ab EA154 d d a 1:2,000.



I c ca a a a a - b dd d Ra B a T
14-3-3 EA154 Rabb Ab d d a 1:500.