

MEK1 Rabbit anti-human MEK1

Catalog NO.: EA219
For research use only.

Overview

Product name MEK1 Rabbit anti-human MEK1

Source Rabbit

Applications WB, IHC

Species reactivity Human, Rat, Mouse

Recommended dilution WesternBlot:1/1000-2000
Immunohistochemistry:1/100-200
NOTE: Optimal dilutions should be determined by the end user.

Immunogen Recombinant Protein

Species Human

Storage PBS + 0.02% sodium azide + 50% H7.4.
Store at -20 °C. Avoid repeated freeze/thaw cycles.

Isotype IgG

Clona-

