



ELK Biotechnology


Kv10.2 Rabbit pAb

Catalog NO.: EA312


For research use only.

Overview

Product name	Kv10.2 Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species reactivity	Human, Rat, Mouse
Recommended dilutions	WesternBlot:1/1000-2000 Immunohistochemistry:1/100-200 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Synthetic peptide
Species	Human
Storage	PBS + 0.02% sodium azide + 50% glycerol at -20°C. Aliquots - avoid repeated freeze-thaw cycles.
Isotype	IG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	60kDa
GeneID Human	238271
Human Swiss-Prot No.	Q8NCM2
Cellular localization	Membrane
Alternative Names	EAG2, H-EAG2, KCNH5
Background	N/A



W b a a 1) Ra Ba T -L M c a P Ma ,
 2)Ra Ba T -H M c a P Ma , 3) M Ba T -
 L M c a P Ma , 4) M Ba T -H M c a
 P Ma , 5) H G2-L M c a P



I c ca a a a a - b dd d Ra Ba T
 K 10.2 (EA312 Rabb Ab d d a 1:200.

I c ca a a a a - b dd d M Ba T
 K 10.2 (EA312 Rabb Ab d d a 1:200.